

FlowChip
MICROSPRAY COLUMN



THE FUTURE OF LC-MS

FEATURES

- High-performance separations at microflow rates
- Robust microflow solution based on proven PicoFrit® technology
- Microspray compatible flow rates
- Easy to use, reproducible, robust universal connectivity
- Portable - Compatible with multiple platforms

HIGH-PERFORMANCE REDEFINED

- Powered by PicoFrit® Technology
- Integral Frit for UHPLC
- Robust, TaperTip™ based
- Inert
- Clog Free
- Sized for Microspray
- Only available from NOI

SPECIFICATIONS

Column ID	150 µm
Bed Length	10 cm
Chemistries	Reposil-PUR (Dr. Maisch)
Flow Rates	1 - 10 µL/min.

As delineated at CPSA Analytics 2019 during the Imagine Summit
www.cpsa-analytics.com/2019/